



'Pivotell Pager' Automatic Pill Dispenser

POWERED BY



Introduction

The new 'Pivotell Pager' Automatic Pill Dispenser offers substantial benefits for those on a regular pill regime who are deaf or have hearing impairment and also have difficulty remembering to take their medication. The dispenser provides important additional support the user, their families and carers.

The 'Pivotell Pager' dispenser has an integral transmitter which triggers an alarm (vibration, sound or flashing light) on any of the Bellman® Visit 868 System receivers to alert the user that it is time to take their medication.

Who is it for?

The Pivotell Pager is for those who are deaf or with hearing impairment who are taking time sensitive medication, who have some confusion (e.g. early dementia), or perhaps learning difficulties and cannot hear the alarm of a Pivotell Mk 3/11 or Pivotell Plus dispenser.

How does it work?

The 'Pivotell Pager' dispenser can hold up to 28 days' supply of medication depending on how frequently medication is taken each day.

When it is time to take medication the Pivotell Pager's internal tray rotates, the alarm sounds and the light flashes. The integral transmitter will activate Bellman® Visit 868 receivers in range which will vibrate or flash to alert the user that it is time to take their medication. The user then simply picks up the dispenser and tilts it to take the medication in the hand or suitable container. This cancels the alarm and flashing light. The vibrating alarm on the receiver can be cancelled by pressing the button. The dispenser will then wait until the next alarm time and repeat the process.

Main features:

- Includes a Bellman® transmitter which alerts Bellman® Visit 868 receivers including:

Bellman® Visit 868
Pager Receiver BE 1470



Bellman® Visit 868
Wrist Receiver B8101



Bellman® Visit 868
Flash Receiver BE1440



- Daily Alarm Times - 1-28
- Compartments - 29 (for easy pharmacy filling one compartment is left empty)
- Two or more compartments can be programmed to alarm at the same time (for those with many or large pills)
- Start delay function - to allow pharmacists or family members to easily fill in advance away from the user's home without having to take out the batteries or set the Doses to '0' (to prevent it from alarming and rotating) and then re-set immediately before delivery to the user.
- Automatic Daylight Saving adjustment (GMT / BST)
- Pharmacy label carrier and label case supplied (to comply with Royal Pharmacy Society of Great Britain labelling requirements)
- Day / time discs supplied (for 1 / 2 / 3 / 4 times daily dosage and one blank disc for five or more doses per day.
- Battery life - 12 months
- Battery level meter and warning indicator
- Secure battery case cover to deter tampering (screw provided)
- Spare trays, lids and time discs available for off-site filling by family members or pharmacists.
- Disposable tray system also available for off-site filling by pharmacists.
- A PivoteLL Pill Dispenser Tipper is available for those with limited strength or wrist dexterity and those who might forget to turn the dispenser upright after dispensing their medication.



- Four different alarm sounds
- Flashing light for those with hearing impairment
- Lockable - only the medication due can be accessed by the user
- The medication tray rotation is disabled if the dispenser is left upside down to avoid pill spillage.

Dimensions:

Diameter – 190mm (7 ½ ins)

Hight 56mm (2 ½ ins)

Weight: 489 gm

PivoteLL Ltd, P.O.Box 108, Saffron Walden, Essex CB11 4WX Tel: 01799 550979

e-mail: office@pivotell.co.uk www.pivotell.co.uk

Copyright: PivoteLL Ltd